

Diversity. Precision.
Reliability – excellence
across every specialty.

Surgiplex. The Difference You Hold





15 cm / 6"
Ø 11 mm

with air regulation and pilot flame, nickel-plated
for gas pressure up to 500 mm water gauge

SP33-100-32

F - for propane or butane gas

SP33-101-42

N - for natural gas



14 cm / 5 1/2"
Ø 11 mm

with lateral pilot-flame tube, lever switchable
either main flame plus pilot flame burning or pilot flame only

SP33-102-32

F - for propane or butane

Double lever stopcock burner

Wipphahn-Brenner
Mecheros con llave de palanca doble
Becs pour la laboratoire
Becco del gas per laboratorio

Miscellaneous
Verschiedenes
Diversos artículos
Articles divers
Articoli diversi

Micro-Standard

SP33-103-37

8,5 cm / 3 1/4"


Ø 6 mm

with stabilizing chimney, gas regulation with conical valve, flame adjustable, for gas pressure up to 500 mm water gauge, chromium-plated



SP33-103-37

F - for propane or butane gas

				
	SP33-104-03	SP33-105-03	SP33-106-07	SP33-108-08
a x b x c	250 x 70 x 3 mm	425 x 70 x 3 mm	285 x 75 x 7 mm	270 x 65 x 8 mm

					
SP33-108-12	SP33-108-13	SP33-108-15	SP33-109-20	SP33-109-25	SP33-109-30
270 x 65 x 12 mm	230 x 60 x 13 mm	205 x 60 x 15 mm	215 x 65 x 20 mm	215 x 65 x 25 mm	215 x 65 x 30 mm

Cleaning brushes

Reinigungsbürsten
Cepillos de limpieza
Brosses à lavement
Spazzole per pulire

Miscellaneous

Verschiedenes
Diversos artículos
Articles divers
Articoli diversi



SP33-110-00	SP33-110-13	= 215 x 60 x 13 mm
set of 5 brushes	SP33-110-15	= 215 x 65 x 15 mm
SP33-110-13 – SP33-110-28		



SP33-111-30	SP33-111-50	SP33-112-95
300 x 120 x 30 mm	320 x 130 x 50 mm	210 x 95 mm



SP33-113-15

Ø = 325 mm, h= 305 mm
Capacity



SP33-114-15

Lid



Martin

SP33-116-58

58 cm / 22 ⁶/₈"



Sauerbruch

SP33-117-55

55 cm / 21 ⁴/₈"

15 l



SP33-113-15